

RNeasy[®] Mini Kit, Part 2

The RNeasy Mini Kit (cat. nos. 74104 and 74106) can be stored at room temperature (15–25°C) for at least 9 months.

For more information, additional and more detailed protocols, and safety information, please refer to the *RNeasy Mini Handbook*, which can be found at www.qiagen.com/handbooks.

For technical assistance, please call toll-free 00800-22-44-6000, or find regional phone numbers at www.qiagen.com/contact.

On-column DNase digestion

Notes before starting

- If using the RNase-Free DNase Set for the first time, prepare DNase I stock solution by injecting 550 μ l RNase-free water into the DNase I vial using an RNase-free needle and syringe. Mix gently by inverting the vial. Do not vortex.
 - For long-term storage of DNase I stock solution, divide it into single-use aliquots and store at –20°C for up to 9 months. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze aliquots after thawing.
1. Add 350 μ l Buffer RW1 to RNeasy column, close lid, centrifuge for 15 s at $\geq 8000 \times g$ ($\geq 10,000$ rpm). Discard flow-through.
 2. Add 10 μ l DNase I stock solution (see above) to 70 μ l Buffer RDD. Mix by gently inverting the tube. Centrifuge briefly.
 3. Add DNase I incubation mix (80 μ l) directly to RNeasy column membrane, and place on benchtop (20–30°C) for 15 min.
 4. Add 350 μ l Buffer RW1 to RNeasy column, close lid, centrifuge for 15 s at $\geq 8000 \times g$. Discard flow-through. Continue with step 5 of “RNA purification from cells/tissue samples” in *Quick-StartProtocol RNeasy Mini Kit, Part 1*, or step 4 of “RNA cleanup” (below).

January 2011



RNA cleanup

Notes before starting

- Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
- 1. Adjust the sample to a volume of 100 μ l with RNase-free water. Alternatively, follow steps in “DNase digestion of RNA before RNA cleanup” in Appendix E of *RNeasy Mini Handbook*. Add 350 μ l Buffer RLT, and mix well.
- 2. Add 250 μ l ethanol (96–100%) to the diluted RNA, and mix well by pipetting. Do not centrifuge. Proceed immediately to step 3.
- 3. Transfer the sample (700 μ l) to an RNeasy Mini spin column placed in a 2 ml collection tube (supplied). Close the lid. Centrifuge for 15 s at $\geq 8000 \times g$. Discard the flow-through.

Optional: If performing optional on-column DNase digestion, follow steps 1–4 of “On-column DNase digestion” (above) after this step.

- 4. Add 500 μ l Buffer RPE to the RNeasy spin column. Close the lid. Centrifuge for 15 s at $\geq 8000 \times g$ to wash the membrane. Discard the flow-through.
- 5. Add 500 μ l Buffer RPE to the RNeasy spin column. Close the lid. Centrifuge for 2 min at $\geq 8000 \times g$ to wash the membrane.

Optional: Place the RNeasy spin column in a new 2 ml collection tube (supplied). Close the lid, and centrifuge at full speed for 1 min.

- 6. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50 μ l RNase-free water directly to the spin column membrane. Close the lid, and centrifuge for 1 min at $\geq 8000 \times g$ to elute the RNA.
- 7. If the expected RNA yield is $>30 \mu$ g, repeat step 6 using another 30–50 μ l of RNase-free water. Alternatively, use the eluate from step 6 (if high RNA concentration is required). Reuse the collection tube from step 6.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, RNeasy® (QIAGEN Group).
1067548 01/2011
© 2011 QIAGEN, all rights reserved.

