

**REED LABORATORY – CONFOCAL IMMUNOFLUORESCENCE**  
**UIC College of Dentistry**

- [0] Fix tissue in 4% PFA overnight  
Embed in Paraffin and section
- [1] **DEPARANNIZING – TIME: 1HR**  
Fill each of the containers  
Xylenes – 15 min  
100% ETOH – 5 min  
100% ETOH – 5 min  
95% ETOH – 5 min  
70% ETOH – 5 min  
1x PBS – 15 min
- [2] **SODIUM BOROHYDRIDE – TIME: 2 X 30 MIN**  
Make fresh, it will start to bubble:  
1x PBS – 15ml  
NaBH<sub>4</sub> – 0.075g (sigma 452882 – 100g)  
Place slides in mailer  
Change after 30 min
- [3] **CITRATE BUFFER – TIME: 10 MIN HEATED, 30 MIN ROOM TEMP**  
Preheat Heat tap water to a boil for 6 min in microwave  
Place on hotplate and bring to boil  
Fill mailer with 0.01 M Citrate Buffer  
Place in boiling water, and take off of the heat source  
Take out of heated water after 10 min  
Allow to cool to room temp for 30 min  
Citrate Buffer, Ph 6 STOCK SOL:  
1.05g citric acid, monohydrate (Fisher BP339,500)  
500ml Dh<sub>2</sub>O  
ADJUST TO Ph 6 WITH 5m NaOH (5M → 20g NaOH in 100 ml Dh<sub>2</sub>O)
- [4] **TRITON WASH – TIME: 3 X 5 MIN**  
Wash in the mailer  
0.5% (v/v) Triton-x100  
To make this, cut the pipette tip off  
Slightly warm PBS before adding Triton  
Eject pipette tip right into PBS
- [5] **MEOH WASH -TIME: 3 MIN**  
Fill mailer with -20 MEOH  
Place slides into MEOH
- [6] **SERUM BLOCKING – TIME: 1 HR**  
Fill mailer with blocking solution, transfer slide right from MEOH into blocking solution  
5% Serum Blocking Solution  
1x PBS – 19 ml  
Animal Serum – 1 ml

We almost exclusively use DONKEY SERUM (Sigma, pn)  
Save 0.5 ml of serum for the primary antibody step

- [7] **PRIMARY ANTIBODY – TIME: O/N**  
Transfer slides to hydrator  
Wipe off excess solution around tissue  
Pipette 250 uL of antibody solution onto slide  
Cover each slide with a slip of PARAFILM  
Cover hydrator, and place in 4 degree refrigerator overnight  
    Primary Antibody Solution [PER SLIDE]  
        5% Serum Blocking Solution - 6.25 µl  
        1x PBS -191.25 µl  
        0.5% Triton – 50 µl  
        Antibody1 (1:200 dilution) – 2.5 µl  
        Antibody2 (1:200 dilution) – 2.5 µl

## DAY 2

- [8] **TRITON WASH – TIME: 3 X 5MIN**
- [9] **SECONDARY ANTIBODY STEP – TIME: 2 HRS**  
*Turn off the overhead lights*  
Transfer slides to hydrator  
Wipe off excess solution around tissue  
Pipette 250 µl of secondary antibody solution onto slide  
Cover hydrator, and place in drawer out of light  
    Secondary Antibody Solution [PER SLIDE]  
        Animal Serum - 5 µl  
        1x PBS - 243.75 µl  
        Secondary Antibody1 (1:200 dilution) – 1.25 µl  
        Secondary Antibody2 (1:200 dilution) – 1.25 µl
- [10] **PBS WASH – TIME: 3 X 5 MIN**
- [11] **DAPI STAIN – TIME: 5 MIN**  
Transfer slides to hydrator  
Wipe off excess solution around tissue  
Pipette 250 uL of DAPI solution onto slide  
    DAPI Solution  
        1x PBS – 250 µl  
        DAPI – 0.25 µl
- [12] **COVER SLIPPING**  
Cover with prolong gold antifade (p36930, Invitrogen)  
Pipette 100 µl on to the coverslip  
Cut the tip of the pipette

STORE THE SLIDES IN THE 4 DEGREE FRIDGE O/N  
SEAL WITH NAIL POLISH THE NEXT DAY